CORONEX-COVID 19 Ver 2.0.

Multiplex RT-qPCR Detection Kit [1000 rxn]-

Ref:DS-CSC19-M1000

DESCRIPTION

CORONEX-COVID 19 (Ver 2.0) qPCR kit provides detection in a single tube, including two independent genes of 2019 Sars-CoV2 and an internal control which targets the human RNaseP (RNP) gene to assess specimen quality and reaction efficiency. Specific primers and probes were designed for the detection of conserved region of 2019 nCoV's 'Orflab' and 'N' gene, respectively, avoiding nonspecific interference of SARS2003 and BatSARS-like virus strains. Internal control (RNaseP gene) provides a nucleic acid extraction procedural control and a secondary negative control.

PRODUCT CONTENTS

Components	Ingredient	Cap number	Amount
CORONEX-Covid19 DS Mix E	RT-qPCR Master mix	1	10X1250 µl
CORONEX-Covid19 DS PP1	Orflab, N and RNP gene primer probe mix	2	5X 500 µl
CORONEX-Covid19 DS PTC*	Positive control for Orflab, N and RNP genes	3	200 µl
No-Template Control (NTC)	No-Template Control	4	200 µl

*tup no 1 and 2 contents are prepared more than 10%.

The "CORONEX-Covid19 DS PTC" contained in the kit is a synthetic DNA prepared to be used as the positive control. It contains the target regions of all genes (Orflab, N and RNP) that are amplified with the primer-probe sets included within the kit. It is not isolated from the virus; it does not have the entire virus genome and has no infectious properties.

STORAGE AND SHELF LIFE

All reagents should be stored at -25 °C to -15 °C with protection from light. The reagents are stable for 14 months when stored at the recommended condition. The expiration date will not change if the kit is opened and stored at the recommended condition. The expiration date will not change if the kit is transported with ice-packs for 4 days and/or treated with 6 freeze-thaw cycles.

INSTRUMENTS

Real time PCR instrument with FAM and HEX channels, such as Magnetic Induction Cycler BMS Mic qPCR, Qiagen Rotor-Gene Q, Biorad CFX96, Montania 4896 Real-Time PCR and LightCycler 480 (Roche).

SAMPLING

Suitable specimen type: upper respiratory specimen (including nasal swabs, nasopharyngeal swabs/aspirates/washes and sputum) and lower respiratory specimen (including respiratory aspirates, bronchial washes and bronchoalveolar lavage fluids).

PROTOCOL

1. Specimen Preparation

The samples should be extracted according to the manufacturer's protocol of viral RNA extraction kits. The extracted RNA can be directly used as nucleic acid template for the reaction.

2. Reagent Preparation

Thaw the required reagents, mix by shaking, and centrifuge briefly before use. Prepare the mixture in a RNase-free centrifuge tube as follows [1]:

NOTE: All work should be carried out on ice/cold plate.

Components	Cap number	Volume (µl Per Reaction)
CORONEX-Covid19 DS Mix E	1	12.5 µl
CORONEX-Covid19 DS PP1	2	2.5 µl
Total Volume	-	15 µl

[1] Calculate the number of reaction tubes (sample number + positive control + negative control). It is recommended to set both negative and positive controls for each test.

Mix the above mixture thoroughly, and make aliquots of 15 µl into different PCR reaction tubes.

3. Template Addition

Add 5 µl of extracted RNA from specimen to different PCR reaction tubes which contained 15 µl of PCR mix. The total volume is 20 µl. Cap the reaction tubes tightly, centrifuge them at low speed.

CORONEX-COVID 19 Ver 2.0.

4. qRT-PCR Amplification

Put the reaction tubes on a PCR instrument, setup and run the following cycling protocol:

Step	Temperature	Time	Cycles
Reverse Transcription	48 °C	20 min	1
Pre-denaturation	95 °C	2 min	1
PCR Cycles	95 °C	5 sec	35
	60 °C (Read)	10 sec	

Settings of detection fluorescence: Orflab and N genes (FAM) and RNP (HEX). Please set the internal reference parameter of fluorescence of the instrument to "None". For example: for ABI series instruments, please set "Passive Reference" to "None".

5. Data Analysis

For Biorad real time PCR instruments;

FAM and HEX:

Threshold should be set to: auto threshold

Analysis Profile for BMS MIC instrument;

Set analysis parameters as follows;

FAM and HEX:

Method: Dynamic

Ignore Cycles Before :5

Exclusion: Extensive,

Fluorescence Cut Off Level: %15 should be on.

Analysis Profile for RotorGeneQ instrument;

When starting a run, make sure Green and Yellow channels are on the channel setup list. Select Gain Optimization and check box 'Perform Optimization Before 1st Acquisition' for both green and yellow channels.

FAM and HEX:

Select dynamic tube.

Ignore first 10 cycles.

Outlier Removal should be set to '5%'.

Eliminate cycles before should be set to:0

Threshold should be set to: 0.03

INTERPRETING PATIENTS' TEST RESULTS

If the criteria of QUALITY CONTROL is met, analyze the data of sample as follows:

- Please follow sigmoidal curves in wells which you observe Ct values.

- The cut-off value in the FAM channel should be evaluated according to the table below#.

In Peltier Block Devices (Biorad CFX, Montania etc)	≤ 35
In Rotor System devices (Qiagen Rotorgene, BMS Mic etc)	≤ 33

#In CDC statements and internationally published literature, it has been reported that passive virus residues can be detected in advanced cycles (1-7).

- If the Ct value of HEX channel is ≤ 33, analyze the results according to the following table:

Orflab and N Gene (FAM)	RNP Gene (HEX)	Interpretation
+	+	The patient is SARS-CoV-2 positive.
-	+	The patient is SARS-CoV-2 negative.
-	-	The result is uncertain. The test should be repeated.*1
+	-	The result is uncertain. The test should be repeated.*1

*1. If the results remain same, repeat the specimen collection/RNA extraction and perform the test again. It is kindly recommended to repeat specimen collection after two days.

References

- Jefferson T., Spencer E., Brassey J., Heneghan C. Viral cultures for COVID-19 infectivity assessment. medRxiv, 2020, doi: <https://doi.org/10.1101/2020.08.04.20167932>
- Bullard J, Durst K, Funk D, Strong JE, Alexander D, Garnett L et al. Predicting Infectious SARS-CoV-2 From Diagnostic Samples. Clinical Infectious Diseases, 2020, doi: 10.1093/cid/ciaa638.
- La Scola B, Le Bideau M, Andreani J, et al. Viral RNA load as determined by cell culture as a management tool for discharge of SARS-CoV-2 patients from infectious disease wards. Eur J Clin Microbiol Infect Dis. 2020;39(6):1059-1061. doi:10.1007/s10096-020-03913-9
- Yu F, Yan L, Wang N, Yang S, Wang L, Tang Y, Gao G, Wang S, Ma C, Xie R, Wang F, Tan C, Zhu L, Guo Y, Zhang F. Quantitative Detection and Viral Load Analysis of SARS-CoV-2 in Infected Patients, Clinical Infectious Diseases, 2020; 71(15): 793–798, <https://doi.org/10.1093/cid/ciaa345>
- Mina M (Harvard T.H. Chan School of Public Health). Your Coronavirus Test Is Positive. Maybe It Shouldn't Be. The New York Times, 2020. <https://www.nytimes.com/2020/08/29/health/coronavirustest.html>
- U.S.A. Centers for Disease Control and Prevention (CDC) report: Duration of Isolation and Precautions for Adults with COVID-19. 2020. <https://www.cdc.gov/coronavirustest/2019-ncov/hcp/duration-on-solat-on.html>
- Drew RJ, O'Donnell S, LeBlanc D, McMahon M, Natin D. The importance of cycle threshold values in interpreting molecular tests for SARS-CoV-2 [published online ahead of print, 2020 Jul 10]. Diagn Microbiol Infect Dis. 2020;98(3):115130. doi:10.1016/j.diagmicrobio.2020.115130

Contact Us:

DS Bio and Nano Technologies Product Tracing and Tracking Corporation.
Address: İçişleri Bakanlığı Mah. Muhsin Yazıcıoğlu Cad. 45/7 Ankara 06530 Turkey

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

Declaration of Conformity

In accordance with EC Directing 98/79/EC relating to in vitro diagnostic medical device(test kit). We herewith declare that the under-mentioned test kit, in view of its design and type of construction, meets the requirements of the above EC Directive.

If device(test kit) is modified without the agreement of the under-signed, this declaration becomes invalid.

Description of Device : CORONEX-Covid19 Test Kit
Relevant EC Directive : In Vitro Diagnostic Directive (IVDD) 98/79/EC (Annex III)
Classification : IVD Other
Harmonized Standarts :ISO/IEC 17025:2005, IVDD 98/79/EC (Annex III), OIE Guideline
Manufacturer : DS Bio and Nanotechnology Product Tracing and Tracking Co
Isçi blokları Mah. Muhsin Yazıcıoglu Cad. 45/7
Ankara 06530 Turkey

Authorized Representative :Dinçer SEN

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has been found in Compliance with requirements of
Environmental Management System

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(Scope in English: Nucleotide Production, Oligo Production, DNA-RNA Production,
DNA-RNA Design, Bioinformatics, Software Business, Software Development.)

Certificate No. : EMS/08134/0918
Original Certificate Date : 24 - September - 2018
Issue Date : 24 - September - 2018
Expiry Date : 23 - September - 2021

To check this certificate status visit:
"<http://uasl.uk.com/certifiedorganization.html>"

Authorised Signature

Quality Control Certification

UK Office: 1929, Chynoweth House,
Trevisson Park, Truro-TR48UN, Cornwall, UK
India Office: 2nd Floor, Aman Market,
Narela Mandi, Delhi - 110 040, India

"Quality Control Certification (QCC)" accredited by "UASL, England, UK". This certificate remains the property of "QCC" to whom it must be returned on request.

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is in Compliance with
Testing and Calibration of Laboratories

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(Scope in English: Nucleotide Production, Oligo Production, DNA-RNA Production,
DNA-RNA Design, Bioinformatics, Software Business, Software Development.)

Certificate No. : TCL/08086/0918
Original Certificate Date : 24 - September - 2018
Issue Date : 24 - September - 2018
Expiry Date : 23 - September - 2021

To check this certificate status visit:
"http://qccertification.com/Client.aspx"

Authorised Signature

For Quality Control Certification

Head Office: 2nd Floor, Aman Market,
Narela Mandi, Delhi - 110 040, India

UK Office: 1929, Chynoweth House,
Trevissome Park, Truro - TR4 8UN,
Cornwall, UK

Quality Control Certification accredited by UASL, UK
This certificate doesn't provide the certified organisation with immunity from its legal obligations.
This certificate remains the property of QC Certification to whom it must be returned on request.

APPLICATION PROTOCOL

Program the qPCR devices as follows and add the reagents to the qPCR tubes in the order specified below, close the tubes, place them into the qPCR device and start the run.

Reaction setup			qPCR Program		
Component	Add order	Reaction	Cycle Number	Temperature	Duration
CORONEX-Covid19 DS Mix-E	1	5 µL	1	45 °C	10 min.
CORONEX-Covid19 DS FP1	2	0.2 µL	1	95 °C	2 min.
CORONEX-Covid19 DS RP1	3	0.2 µL	40	95 °C	5 sec.
CORONEX-Covid19 DS P1	4	0.2 µL		60 °C	20 sec.
CORONEX-Covid19 DS C1	7	0.1 µL			
CORONEX-Covid19 DS C2	8	0.2 µL			
CORONEX-Covid19 DS PCT1 (positive control)	6	1 µL			
Template		3.5 µL			
TOTAL REACTION VOLUME		10µL		FAM	

WARNING: Interpretation of Results: 1) Examine the shape of the amplification curves obtained in the FAM channels and record non-sigmoidal curves as negative. 2) Calculate the number of threshold cycles (Cq). 3) Record the result as negative if $38 \leq Cq$ and as positive if $Cq < 38$.

Intended Use and Test Principle: Kit is used for testing 2019-nCoV causing Coronavirus Disease 2019 (COVID-19). The kit is applied to nucleic acid isolates from nasopharyngeal aspirate/lavage, bronchoalveolar lavage, nasopharyngeal swab, oropharyngeal swab and sputum samples. Rapid diagnosis with the kit is achieved both conventional PCR (cPCR) and real-time PCR (qPCR) targeting 2019-nCoV ORF (open reading frame) viral translational region. This region-specific oligonucleotide primer and prob sets gives positive results only with 2019-nCoV. In the 2019-nCoV routine PCR screening, CORONEX-COVID 19 test kit, primer and prob set is applied; if the result is 2019-nCoV ORF (open reading frame) viral translational region positive, the patient result can be interpreted as infection positive. If the result is negative, patient result can be interpreted as disease free.

Analytical Specification: Kit is validated Bio-Rad CFX384 and Biomolecular Systems MIC qPCR machine. The kit's testing limit (LOD) is the lowest viral copy that can be tested with 99% probability.

Warnings: Master stock reagents should be kept on the cold block during cPCR and qPCR set up. If possible, all PCR set up should be performed on the cold block.

Kit components should be mixed by gently shaking before use. The micropipettes used for pipetting cPCR and qPCR mixes, positive control and patient samples should be separated.

Kit content (Shelf Life: 12 months) Storage: -20 °C;
Transport: +2-8 °C, 200 rxn.

Intended Use/Content/Target	Quantity	Consumption per reaction
CORONEX-Covid19 DS Mix-E	1 X 1000 µL	5 µL
CORONEX-Covid19 DS FP1	1 X 80 µL	0.4 µL
CORONEX-Covid19 DS RP1	1 X 80 µL	0.4 µL
CORONEX-Covid19 DS P1	1 X 80 µL	0.4 µL
CORONEX-Covid19 DS C1	1 X 20 µL	0.1 µL
CORONEX-Covid19 DS C2	1 X 40 µL	0.2 µL
CORONEX-Covid19 DS PCT1 (positive control)	1 X 10 µL	1 µL

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is in Compliance with
Medical Laboratories
Requirements for Quality and Competence

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(Scope in English: Nucleotide Production, Oligo Production, DNA-RNA Production,
DNA-RNA Design, Bioinformatics, Software Business, Software Development.)

Certificate No. : MLRQC/08087/0918
Original Certificate Date : 24 - September - 2018
Issue Date : 24 - September - 2018
Expiry Date : 23 - September - 2021

To check this certificate status visit:
"http://qccertification.com/Client.aspx"

Authorised Signature

For Quality Control Certification

Head Office: 2nd Floor, Aman Market,
Narela Mandi, Delhi - 110 040, India

UK Office: 1929, Chynoweth House,
Trevisson Park, Truro - TR4 8UN,
Cornwall, UK

Quality Control Certification accredited by UASL, UK
This certificate doesn't provide the certified organization with immunity from its legal obligations.
This certificate remains the property of QC Certification to whom it must be returned on request.

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is in Compliance with
Risk Management
ISO 31000:2018

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(Scope in English: Nucleotide Production, Oligo Production, DNA-RNA Production,
DNA-RNA Design, Bioinformatics, Software Business, Software Development.)

Certificate No. : RM/08121/0918
Original Certificate Date : 24 - September - 2018
Issue Date : 24 - September - 2018
Expiry Date : 23 - September - 2021

To check this certificate status visit:
"<http://qccertification.com/Client.aspx>"

Authorised Signature

For Quality Control Certification

Head Office: 2nd Floor, Aman Market,
Narela Mandi, Delhi - 110 040, India

UK Office: 1929, Chynoweth House,
Trevissome Park, Truro - TR4 8UN,
Cornwall, UK

Quality Control Certification accredited by UASL, UK
This certificate doesn't provide the certified organisation with immunity from its legal obligations.
This certificate remains the property of QC Certification to whom it must be returned on request.