

THO LISTER OF

1032159

**Each vial contains:** Remdesivir 100 mg. **Dosage:** As directed by the Physician.

See enclosed prescribing information for reconstitution instructions and complete information on dosage and administration.

Store below 30°C until required for use. Reconstitute with 19 mL of sterile water for injection and diluted into 0.9% saline prior to administration by intravenous infusion.

After reconstitution, the storage time before administration should not exceed 4 hours at room temperature (20°C to 25°C) or 24 hours at refrigerated temperature (2°C to 8°C).

### Preservative Free

Note: Parenteral products should be inspected visually for particulate matter prior to administration.

KEEP OUT OF THE REACH AND SIGHT OF

For more information please scan the code >>>>>





# **DESREM**<sup>™</sup>

Lyophilized Powder for Injection for IV Infusion

Discard Unused Portion Single-Dose Vial Sterile

POM Schedule 2 PP

NOT FOR DISTRIBUTION IN US, CANADA OR EU





DESREM is manufactured under a license from Gilead Sciences, Inc.

To be distributed and used in licensed countries only. See packaged insert for licensed countries.

Manufactured by:

Mylan Laboratories Limited,

At: Plot No. 42 to 52, Survey No. 166, 171, 172 & 177, TSIIC, Phase -III, IDA, Pashamylaram (V), Patancheru (M), Sangareddy District, Telangana State – 502307.

Code No.: TS/DRUGS/SGY/2020-65910 TM - Trade Mark under registration



# Remdesivir for Injection 100 mg/vial

# **DESREM**<sup>™</sup>

Lyophilized Powder for Injection for IV Infusion

Discard Unused Portion Single-Dose Vial Sterile

POM Schedule 2 PP

NOT FOR DISTRIBUTION IN US, CANADA OR EU







## ARTWORK DETAIL LABEI

ARI WORK DETAIL LABEL					
Product	DESREM - Remdesivir for Injection 100 mg/vial				
Buyer/Country	Mylan / RoW (General Export) - (GLAND)	Component	Carton		
Dimension	48 x 48 x 88 mm			Pack	1's
New Item Code	1032159	Old Item Code	NA		
Colour Shades	PMS 7718 C PMS 7516 C	PMS 306 C	Black	No. of Colours	4
Change Control No.	NA				
Design/Style	Top Open, Reverse Tuck.				
Substrate	300 GSM CFBB with spot aqueous varnish. Board thickness in microns 465.				
Special Instructions	Printing clarity to be clear & sharp.				
Autocartonator Requirements	NA				
	r: Before processing, please ensure that the ARTV se of any FONTS/DESIGN are Mis-matching with the				

SOP - 000484829 Rev. 0; Ver. 1; Date: 02-10-2020

MAKE ANY CHANGE TO THE ARTWORK WITHOUT WRITTEN INSTRUCTIONS FROM PDC.